

# Package ‘synapter’

December 28, 2024

**Type** Package

**Title** Label-free data analysis pipeline for optimal identification and quantitation

**Version** 2.31.0

**Author** Laurent Gatto, Nick J. Bond, Pavel V. Shliaha and Sebastian Gibb.

**Maintainer** Laurent Gatto <laurent.gatto@uclouvain.be>  
Sebastian Gibb <mail@sebastiangibb.de>

**Depends** R (>= 3.1.0), methods, MSnbase (>= 2.1.2)

**Imports** RColorBrewer, lattice, qvalue, multtest, utils, tools,  
Biobase, Biostrings, cleaver (>= 1.3.3), readr (>= 0.2),  
rmarkdown (>= 1.0)

**Suggests** synapterdata (>= 1.13.2), xtable, testthat (>= 0.8), BRAIN,  
BiocStyle, knitr

**Description** The synapter package provides functionality to reanalyse label-free proteomics data acquired on a Synapt G2 mass spectrometer. One or several runs, possibly processed with additional ion mobility separation to increase identification accuracy can be combined to other quantitation files to maximise identification and quantitation accuracy.

**License** GPL-2

**URL** <https://lgatto.github.io/synapter/>

**VignetteBuilder** knitr

**biocViews** ImmunoOncology, MassSpectrometry, Proteomics, QualityControl

**RoxygenNote** 6.1.0

**git\_url** <https://git.bioconductor.org/packages/synapter>

**git\_branch** devel

**git\_last\_commit** 4502012

**git\_last\_commit\_date** 2024-10-29

**Repository** Bioconductor 3.21

**Date/Publication** 2024-12-27

## Contents

synapter-package . . . . .	2
createUniquePeptideDbRds . . . . .	3
Deprecated . . . . .	4
estimateMasterFdr . . . . .	4
inspectPeptideScores . . . . .	5
makeMaster . . . . .	6
MasterFdrResults-class . . . . .	8
MasterPeptides-class . . . . .	9
plotFragmentMatching . . . . .	10
requantify,MSnSet-method . . . . .	11
rescaleForTop3,MSnSet,MSnSet-method . . . . .	12
Synapter . . . . .	13
synapterGuide . . . . .	26
synapterPlgsAgreement,MSnSet-method . . . . .	27
synapterTinyData . . . . .	28
synergise . . . . .	29
synergise2 . . . . .	32
<b>Index</b>	<b>36</b>

---

synapter-package	<i>Combine label-free data for optimal identification and quantitation</i>
------------------	--

---

## Description

The synapter package provides functionality to reanalyse label-free proteomics data acquired on a Synapt G2 mass spectrometer. One or several runs, possibly processed with additional ion mobility separation to increase identification accuracy can be combined to other quantitation files to maximise identification and quantitation accuracy.

## Details

Two pipelines of varying flexibility are proposed the preform data analysis: (1) the `synergise1` and `synergise2` function are single entry functions for a complete analysis and (2) low level, step-by-step data processing can be achieved as described in `?Synapter`.

A high-level overview of the package and how to operate it can be found in the vignette, accessible with `synapterGuide()`. Detailed information about the data processing can be found in the respective function and class manual pages appropriately referenced in the vignette.

For questions, use the Bioconductor mailing list or contact the author. The vignette has a section with details on where/how to get help.

## Author(s)

Laurent Gatto, Pavel V. Shliha, Nick J. Bond and Sebastian Gibb

Maintainer: Laurent Gatto <lg390@cam.ac.uk> and Sebastian Gibb <mail@sebastiangibb.de>

## References

*Improving qualitative and quantitative performance for MSE-based label free proteomics*, N.J. Bond, P.V. Shliaha, K.S. Lilley and L. Gatto, J Proteome Res. 2013 Jun 7;12(6):2340-53. doi: 10.1021/pr300776t. PubMed PMID: 23510225.

*The Effects of Travelling Wave Ion Mobility Separation on Data Independent Acquisition in Proteomics Studies*, P.V. Shliaha, N.J. Bond, L. Gatto and K.S. Lilley, J Proteome Res. 2013 Jun 7;12(6):2323-39. doi: 10.1021/pr300775k. PMID: 23514362.

---

createUniquePeptideDbRds

*Create an RDS file for the 'Unique Peptides Database'*

---

## Description

This function creates an RDS file to store the “Unique Peptides Database” for future runs of [synergise](#) or [Synapter](#).

## Usage

```
createUniquePeptideDbRds(fastaFile, outputFile = paste0(fastaFile, ".rds"),
  missedCleavages = 0, IisL = FALSE, verbose = interactive())
```

## Arguments

fastaFile	file path of the input fasta file
outputFile	file path of the target RDS file; must have the file extension ".rds"
missedCleavages	Number of maximal allowed missed cleavages. Default is 0.
IisL	If TRUE Isoleucin and Leucin are treated as identical. In this case sequences like "ABCI", "ABCL" are removed because they are not unique. If FALSE (default) "ABCI" and "ABCL" are reported as unique.
verbose	If TRUE a verbose output is provided.

## Author(s)

Sebastian Gibb <mail@sebastiangibb.de>

## See Also

[Synapter](#) for details about the cleavage procedure.

## Examples

```
## Not run:
createUniquePeptideDbRds("uniprot.fasta", "uniprot.fasta.rds")

## End(Not run)
```

---

Deprecated	<i>synapter Deprecated and Defunct</i>
------------	--

---

### Description

The function, class, or data object you have asked for have been deprecated or made defunct. Deprecated functions are provided for compatibility with older versions of the synapter package only, and will be defunct at the next release. Defunct functions have been removed.

Deprecated:

Defunct: synapterGUI; use the replacement [synergise](#)

---

estimateMasterFdr	<i>Computes FDR for all possible final peptide combinations</i>
-------------------	---

---

### Description

This function takes all possible combination of pepfiles of length greater or equal than 2 and computes the number of estimated incorrect peptides, the number of unique peptides, the number of unique proteotypic peptides and the false discovery rate after merging for each combination. The best combination has an fdr lower than masterFdr and the highest number of unique (proteotypic) peptides.

### Usage

```
estimateMasterFdr(pepfiles, fastafile, masterFdr = 0.025, fdr = 0.01,
  proteotypic = TRUE, missedCleavages = 0, IisL = FALSE,
  maxFileComb = length(pepfiles), verbose = interactive())
```

### Arguments

pepfiles	A list of vector of final peptide filenames.
fastafile	A character with the fasta filename.
masterFdr	A numeric indicating the maximum merged false discovery to be allowed.
fdr	Peptide FDR level for individual peptide files filtering.
proteotypic	Logical. Should number proteotypic peptides be used to choose best combination and plot results or total number of unique peptides.
missedCleavages	Number of maximal missed cleavage sites. Default is 0.
IisL	If TRUE Isoleucin and Leucin are treated as identical. In this case sequences like "ABCI", "ABCL" are removed because they are not unique. If FALSE (default) "ABCI" and "ABCL" are reported as unique.
maxFileComb	A numeric to limit the accepted file combinations to reduce computation time. Default is length(pepfiles) meaning no limit set.
verbose	Should progress messages be printed?

### Details

The false discovery rate for the master (merged) file is calculated by summing the number of estimated false discoveries for each individual final peptide file (number of unique peptides in that file multiplied by `fdr`) divided by the total number of unique peptides for that specific combination.

The function returns an instance of the class "[MasterFdrResults](#)".

### Value

An instance of class "[MasterFdrResults](#)". See details above.

### Author(s)

Laurent Gatto

### References

Bond N. J., Shliaha P.V., Lilley K.S. and Gatto L., (2013) J. Prot. Research.

### See Also

The [makeMaster](#) function to combine the peptide data as suggested by `estimateMasterFdr` into one single *master* peptide file.

The vignette, accessible with `synapterGuide()` illustrates a complete pipeline using `estimateMasterFdr` and `makeMaster`.

---

`inspectPeptideScores`    *Inspect peptide scores.*

---

### Description

This function takes a final peptide file and returns information about the *unique* peptide scores and their number in each peptide matchType (PepFrag1 and PepFrag2) by protein dataBaseType (Random and Regular) category. The function is a lightweight version of [getPepNumbers](#) and [plotPepScores](#) (to be used with Synapter instances) for individual files.

### Usage

```
inspectPeptideScores(filename)
```

### Arguments

`filename`            The name of a final peptide file.

### Value

A table of peptide counts in each peptide matchType \* protein dataBasteType category. Also plots the distribution of respective peptide scores.

**Author(s)**

Laurent Gatto

makeMaster

*Merges final peptide files***Description**

This function combines a list of peptide final peptide files into one single *master* file that is obtained by merging the unique peptides from the filtered original peptide files. Additionally it can combine multiple final fragment files into a fragment library.

**Usage**

```
makeMaster(pepfiles, fragmentfiles, fdr = 0.01, method = c("BH",
  "Bonferroni", "qval"), span.rt = 0.05, span.int = 0.05,
  maxDeltaRt = Inf, removeNeutralLoss = TRUE, removePrecursor = TRUE,
  tolerance = 2.5e-05, verbose = interactive())
```

**Arguments**

pepfiles	A character vector of final peptide file names to be merged.
fragmentfiles	A character vector of final fragment file names to be combined into an fragment library. These files should be from the same runs as the final peptide files used in pepfiles.
fdr	A numeric indicating the peptide false discovery rate limit.
method	A character indicating the p-value adjustment to be used. One of BH (default), Bonferroni or qval.
span.rt	A numeric with the loess span parameter value to be used for retention time modelling.
span.int	A numeric with the loess span parameter value to be used for intensity modelling.
maxDeltaRt	A double value that sets a maximum limit for the retention time deviation between master and slave run to be included
removeNeutralLoss	A logical, if TRUE peptides with neutral loss are removed from the fragment library.
removePrecursor	A logical, if TRUE precursor ions are removed from the fragment spectra.
tolerance	A double value that determines the tolerance used to look for the precursor ions.
verbose	A logical indicating if information should be printed out.

## Details

The merging process is as follows:

1. Each individual peptide final peptide file is filtered to retain (i) non-duplicated unique tryptic peptides, (ii) peptides with a false discovery rate  $\leq$  `fdr` and (iii) proteins with a false positive rate  $\leq$  `fpr`.
2. The filtered peptide files are ordered (1) according to their total number of peptides (for example [P1, P2, P3]) and (2) as before with the first item is positioned last ([P2, P3, P1] in the previous example). The peptide data are then combined in pairs in these respective orders. The first one is called the *master* file.
3. For each (master, slave) pair, the slave peptide file retention times are modelled according to the (original) master's retention times and slave peptides, not yet present in the master file are added to the master file.
4. The final *master* datasets, containing their own peptides and the respective slave specific retention time adjusted peptides are returned as a `MasterPeptides` instance.

The resulting `MasterPeptides` instance can be further used for a complete master vs. peptides/Pep3D analysis, as described in [Synapter](#), [synergise](#) or using the GUI ([synapterGUI](#)). To do so, it must be serialised (using the `saveRDS` function) with a `.rds` file extension, to be recognised (and loaded) as a R object.

When several quantitation (or identification) files are combined as a master set to be mapped back against the individual final peptide files, the second master [P2, P3, P1] is used when analysing the peptide data that was first selected in the master generation (P1 above). This is to avoid aligning two identical sets of peptides (those of P1) and thus not being able to generate a valid retention time model. This is detected automatically for the user.

The two master peptides dataframes can be exported to disk as two csv files with `writeMasterPeptides`. The `MasterPeptides` object returned by `makeMaster` can be saved to disk (with `save` or `saveRDS`) and later reloaded (with `load` or `readRDS`) for further analysis.

The fragment library generation works as follows:

1. Each individual final fragment file is imported and only peptides present in the *master* dataset are used.
2. The fragments are combined based on their precursor ions.
3. The intensities of identical fragments (seen in different runs) is summed and divided by the summed precursor intensity (of the same peptide in different runs).
4. Afterwards the intensities are normalized to the average precursor intensity of the different runs.
5. Finally a [MSnExp](#) object is created.

The fragment library dataframe can be exported to disk as csv file with `writeFragmentLibrary`.

## Value

An instance of class "[MasterPeptides](#)".

## Author(s)

Laurent Gatto, Sebastian Gibb

## References

Shliaha P.V., Bond N. J., Lilley K.S. and Gatto L., in prep.

## See Also

See the [Synapter](#) class manual page for detailed information on filtering and modelling and the general algorithm implemented in the synapter package.

The [estimateMasterFdr](#) function allows to control false discovery rate when combining several peptide files while maximising the number of identifications and suggest which combination of peptide files to use.

The vignette, accessible with `synapterGuide()` illustrates a complete pipeline using `estimateMasterFdr` and `makeMaster`.

---

MasterFdrResults-class

*Class "MasterFdrResults"*

---

## Description

This class stored the results of the

## Objects from the Class

Objects are created with the [estimateMasterFdr](#) function.

## Slots

**all**: Object of class "data.frame" with results of all possible combinations.

**files**: Object of class "character" with file names used in combinations.

**best**: Object of class "data.frame" with results of the best combination.

**masterFdr**: Object of class "numeric" storing the best combination.

## Methods

**bestComb** signature(object = "MasterFdrResults"): ...

**allComb** signature(object = "MasterFdrResults"): ...

**fileNames** signature(object = "MasterFdrResults"): ...

**masterFdr** signature(object = "MasterFdrResults"): ...

**plot** signature(x = "MasterFdrResults", y = "missing"): ...

**show** signature(object = "MasterFdrResults"): ...

## Author(s)

Laurent Gatto <lg390@cam.ac.uk>

## References

Improving qualitative and quantitative performance for MSE-based label free proteomics, N.J. Bond, P.V. Shliaha, K.S. Lilley and L. Gatto, Journal of Proteome Research, 2013, in press.

The Effects of Travelling Wave Ion Mobility Separation on Data Independent Acquisition in Proteomics Studies, P.V. Shliaha, N.J. Bond, L. Gatto and K.S. Lilley, Journal of Proteome Research, 2013, in press.

## Examples

```
## Not run:
library(synapterdata)
loadMaster()
class(master)
master

## End(Not run)
```

---

MasterPeptides-class    *Class "MasterPeptides"*

---

## Description

A class to store the results of [makeMaster](#). This class stored the 2 versions (orders) of the *master* final peptide data.

## Objects from the Class

Objects can be created by calls of the form [makeMaster](#).

## Slots

**masters:** Object of class "list" storing the 2 master data.frame objects.

**pepfiles:** Object of class "character" with the list of final peptide input files.

**fdr:** Object of class "numeric" with the peptide false discovery applied when creating the filter.

**method:** Object of class "character" with the peptide p-value adjustment method. One of BH (default), qval or Bonferroni.

**orders:** Object of class "list" with the numeric vectors specifying the order of pepfiles used to generate the respective masters data.frames.

**fragmentfiles:** Object of class "character" with the list of final fragment input files.

**fragments:** Object of class data.frame storing the combined final fragment data.

**fragmentlibrary:** Object of class "MSnExp" storing the fragment library.

## Methods

**show** signature(object = "MasterPeptides"): to print a textual representation of the instance.

**Author(s)**

Laurent Gatto

**References**

Improving qualitative and quantitative performance for MSE-based label free proteomics, N.J. Bond, P.V. Shliaha, K.S. Lilley and L. Gatto, Journal of Proteome Research, 2013, in press.

The Effects of Travelling Wave Ion Mobility Separation on Data Independent Acquisition in Proteomics Studies, P.V. Shliaha, N.J. Bond, L. Gatto and K.S. Lilley, Journal of Proteome Research, 2013, in press.

---

plotFragmentMatching    *Plot fragment matching results.*

---

**Description**

This method plots the results of the fragment matching procedure ([fragmentMatching](#)). A single plot contains two panels. The upper panel shows the identification fragments and the lower one the MS2 spectrum of the quantitation run. Common peaks are drawn in a slightly darker colour and with a point on the top.

**Arguments**

object	Object of class " <a href="#">Synapter</a> ".
key	character, value to look for.
column	character, name of the column in which plotFragmentMatching looks for key.
...	Further arguments passed to internal functions.

**Methods**

signature(object = "Synapter", key = "character", column = "character", verbose = "logical", ...)

Plots identification fragments against quantitation spectra. The ... arguments are passed to the internal functions. Currently legend.cex, fragment.cex, and most of the [plot.default](#) arguments (like xlim, ylim, main, xlab, ylab, ...) are supported. legend.cex and fragment.cex control the size of the legend and fragments labels (default: 0.5). Please see [par](#) for details about cex. If verbose = TRUE a progress bar is shown.

**Author(s)**

Sebastian Gibb <mail@sebastiangibb.de>

**See Also**

[fragmentMatching](#)

---

requantify, MSnSet-method

*Intensity requantification*


---

## Description

This method tries to remove saturation effects from the intensity counts.

## Usage

```
## S4 method for signature 'MSnSet'
requantify(object, saturationThreshold,
method=c("sum", "reference", "th.mean", "th.median",
"th.weighted.mean"), ...)
```

## Arguments

object	An <a href="#">MSnSet</a> object.
saturationThreshold	double, intensity of an ion (isotope of a given charge state) at which saturation is starting to occur.
method	character, requantification method, please see details section.
...	further arguments passed to internal functions. Currently only <code>onlyCommonIsotopes</code> for <code>method="sum"</code> and <code>requantifyAll</code> for <code>method=c("th.mean", "th.median", "th.weighted.mean")</code> are supported.

## Details

Currently `requantify` supports 3 (5) different requantification methods.

"sum" is the simplest requantification method. All ions of a peptide below the saturation threshold are summed to get the new intensity. This method accept an additional argument, namely `onlyCommonIsotopes`. If `onlyCommonIsotopes=TRUE` (default) all ions that are not seen in all runs are removed and only the common seen ions are summed. In contrast `onlyCommonIsotopes=FALSE` sums all ions regardless they are present in all runs.

In "reference" the run that has the most unsaturated ions in common with all the other runs. If there are more than one run, the most intense is used as reference. The other runs are corrected as follows:

- Find common ions between current and the reference run.
- Divide the intensities of the common ions and calculate the mean of these quotients as a run specific scaling factor.
- Multiply the unsaturated ions of the current run by the scaling factor and replace the saturated ones by the product of the scaling factor and the intensities of their corresponding ions in the reference run.

- Sum the rescaled ion intensities.

The "th.\*" methods are nearly identical. All of them calculate the theoretical isotopic distribution for the given sequence of the peptide. Subsequently the unsaturated ions are divided by their theoretical proportion and the mean/median/weighted.mean (proportions are used as weights) of these intensities are calculated per charge state. The sum of the charge state values is used as requantified intensity for this peptide.

If requantifyAll=FALSE (default) just peptides with at least one saturated ion are requantified (unsaturated peptides are unaffected). If requantify=TRUE all peptides even these where all ions are below saturationThreshold are requantified by their theoretical distribution.

### Value

MSnSet where the assayData are requantified.

### Author(s)

Sebastian Gibb <mail@sebastiangibb.de> and Pavel Shliaha

### References

See discussion on github: <https://github.com/lgatto/synapter/issues/39>

### See Also

MSnSet documentation: [MSnSet](#)

---

rescaleForTop3,MSnSet,MSnSet-method  
*Rescale for TOP3*

---

### Description

This method rescales the intensity values of an [MSnSet](#) object to be suitable for TOP3 quantification.

### Usage

```
## S4 method for signature 'MSnSet,MSnSet'  
rescaleForTop3(before, after, saturationThreshold,  
onlyForSaturatedRuns=TRUE, ...)
```

## Arguments

before	An <a href="#">MSnSet</a> object before requantification.
after	The same <a href="#">MSnSet</a> object as before but after requantification.
saturationThreshold	double, intensity of an ion (isotope of a given charge state) at which saturation is starting to occur.
onlyForSaturatedRuns	logical, rescale just runs where at least one isotope is affected by saturation.
...	further arguments passed to internal functions. Currently ignored.

## Details

If an [MSnSet](#) object was requantified using the method="sum" requantification method (see [requantify,MSnSet-method](#)) TOP3 is not valid anymore because the most abundant proteins are penalised by removing high intensity isotopes.

To overcome this `rescaleForTop3` takes the proportion `isotope/sum(isotopes)` for each requantified peptide and calculates a correction factor by comparing these proportions against the unsaturated isotopes before requantification. The new rescale intensity values of the isotopes are the mean correction factor multiplied with the corrected intensity values (see <https://github.com/lgatto/synapter/issues/39#issuecomment-207987278> for the complete explanation/discussion).

## Value

[MSnSet](#) where the assayData are requantified.

## Author(s)

Sebastian Gibb <[mail@sebastiangibb.de](mailto:mail@sebastiangibb.de)> and Pavel Shliaha

## References

See discussion on github: <https://github.com/lgatto/synapter/issues/39>

## See Also

MSnSet documentation: [MSnSet](#) <https://github.com/lgatto/synapter/issues/39#issuecomment-207987278>

## Description

A reference class to store, manage and process Synapt G2 data to combine identification and quantitation results.

The data, intermediate and final results are stored together in such a ad-hoc container called a class. In the frame of the analysis of a set of 3 or 5 data files, namely as identification peptide, a quantitation peptide and a quantitation Pep3D, and identification fragments and quantitation spectra, such a container is created and populated, updated according to the user's instructions and used to display and export results.

The functionality of the synapter package implemented in the Synapter class is described in the *Details* section below. Documentation for the individual methods is provided in the *Methods* section. Finally, a complete example of an analysis is provided in the *Examples* section, at the end of this document.

See also papers by Shliaha et al. for details about ion mobility separation and the manuscript describing the synapter methodology.

## Usage

```
Synapter(filenamees, master) ## creates an instance of class 'Synapter'
```

## Arguments

filenamees	A named list of file names to be load. The names must be 'identpeptide', 'quantpeptide', 'quantpep3d' and 'fasta' (could be an RDS file created by <code>link{createUniquePeptideDb}</code> and (optional 'identfragments' and 'quantspectra'). identpeptide can be a csv final peptide file (from PLGS) or a saved " <a href="#">MasterPeptides</a> " data object as created by <a href="#">makeMaster</a> if working with <i>master</i> peptide data. To serialise the " <a href="#">MasterPeptides</a> " instance, use the <code>saveRDS</code> function, and file extension <code>rds</code> . This master file could contain a fragment library as well. In this case the 'identfragments' argument would be ignored.
master	A logical that defines if the identification file is a <i>master</i> file. See <a href="#">makeMaster</a> for details about this strategy.

## Details

A Synapter object logs every operation that is applied to it. When displayed with `show` or when the name of the instance is typed at the R console, the original input file names, all operations and resulting the size of the respective data are displayed. This allows the user to trace the effect of respective operations.

**Loading the data:** The construction of the data and analysis container, technically defined as an instance or object of class `Synapter`, is created with the `Synapter` constructor. This function requires 4 or 6 files as input, namely, the identification final peptide csv file, the quantitation final peptide csv file, the quantitation Pep3D csv file (as exported from the PLGS software), the fasta file use for peptide identification, and optional the identification fragments csv file and the quantitation spectra xml file. The fasta file ('fasta') could be an RDS file generated by `link{createUniquePeptideDbRds}`, too. The file names need to be specified as a named list with names 'identpeptide', 'quantpeptide', 'quantpep3d', 'fasta', 'identfragments'

and 'quantspectra' respectively. These files are read and the data is stored in the newly created Synapter instance.

The final peptide files are filtered to retain peptides with matchType corresponding to PepFrag1 and PepFrag2, corresponding to unmodified round 1 and 2 peptide identification. Other types, like NeutralLoss\_NH3, NeutralLoss\_H2O, InSource, MissedCleavage or VarMod are not considered in the rest of the analysis. The quantitation Pep3D data is filtered to retain Function equal to 1 and unique quantitation spectrum ids, i.e. unique entries for multiple charge states or isotopes of an EMRT (exact mass-retention time features).

Then, p-values for Regular peptides are computed based on the Regular and Random database types score distributions, as described in *Käll et al., 2008a*. Only unique peptide sequences are taken into account: in case of duplicated peptides, only one entry is kept. Empirical p-values are adjusted using *Bonferroni* and *Benjamini and Hochberg, 1995* (multtest package) and q-values are computed using the qvalue package (*Storey JD and Tibshirani R., 2003 and Käll et al., 2008b*). Only Regular entries are stored in the resulting data for subsequent analysis.

The data tables can be exported as csv spreadsheets with the writeIdentPeptides and writeQuantPeptides methods.

**Filtering identification and quantitation peptide:** The first step of the analysis aims to match reliable peptide. The final peptide datasets are filtered based on the FDR (BH is default) using the filterQuantPepScore and filterIdentPepScore methods. Several plots are provided to illustrate peptide score densities (from which p-values are estimated, plotPepScores; use getPepNumbers to see how many peptides were available) and q-values (plotFdr).

Peptides matching to multiple proteins in the fasta file (non-unique tryptic identification and quantitation peptides) can be discarded with the filterUniqueDbPeptides method. One can also filter on the peptide length using filterPeptideLength.

Another filtering criterion is mass accuracy. Error tolerance quantiles (in ppm, parts per million) can be visualised with the plotPpmError method. The values can be retrieved with getPpmErrorQs. Filtering is then done separately for identification and quantitation peptide data using filterIdentPpmError and filterQuantPpmError respectively. The previous plotting functions can be used again to visualise the resulting distribution.

Filtering can also be performed at the level of protein false positive rate, as computed by the PLGS application (protein.falsePositiveRate column), which counts the percentage of decoy proteins that have been identified prior to the regular protein of interest. This can be done with the filterIdentProtFpr and filterQuantProtFpr methods. Note that this field is erroneously called a false positive rate in the PLGS software and the associated manuscript; it is a false discovery rate.

**Merging identification and quantitation peptides:** Common and reliable identification and quantitation peptides are then matched based on their sequences and merged using the mergePeptides method.

**Retention time modelling:** Systematic differences between identification features and quantitation features retention times are modelled by fitting a local regression (see the [loess](#) function for details), using the modelRt method. The smoothing parameter, or number of neighbour data points used for local fit, is controlled by the span parameter that can be set in the above method.

The effect of this parameter can be observed with the plotRt method, specifying what = "data" as parameters. The resulting model can then be visualised with the above method specifying what

= "model", specifying up to 3 number of standard deviations to plot. A histogram of retention time differences can be produced with the `plotRtDiffs` method.

To visualise the feature space `plotFeatures` could be used. It generates one or two (if ion mobility is available) plots of retention time vs mass and mass vs ion mobility for each data source, namely, Identification data, Quantitation data and Quantitation Pep3D data.

**Intensity modelling:** Systematic differences between intensities of identification features and quantitation features depending on retention times are modelled by fitting a local regression (see the `loess` function for details), using the `modelIntensity` method. The smoothing parameter, or number of neighbour data points used for local fit, is controlled by the `span` parameter that can be set in the above method.

The effect of this parameter can be observed with the `plotIntensity` method, specifying what = "data" as parameters. The resulting model can then be visualised with the above method specifying what = "model", specifying up to 3 number of standard deviations to plot.

**Grid search to optimise matching tolerances:** Matching of identification peptides and quantitation EMRTs is done within a mass tolerance in parts per million (ppm) and the modelled retention time +/- a certain number of standard deviations. To help in the choice of these two parameters, a grid search over a set of possible values is performed and performance metrics are recorded, to guide in the selection of a 'best' pair of parameters.

The following metrics are computed: (1) the percentage of identification peptides that matched a single quantitation EMRT (called `prcntTotal`), (2) the percentage of identification peptides used in the retention time model that matched the quantitation EMRT corresponding to the correct quantitation peptide in ident/quant pair of the model (called `prcntModel`) and (3) the detailed about the matching of the features used for modelling (accessible with `getGridDetails`) and the corresponding details grid that reports the percentage of correct unique assignments. The detailed grid results specify the number of non matched identification peptides (0), the number of correctly (1) or wrongly (-1) uniquely matched identification peptides, the number of identification peptides that matched 2 or more peptides including (2+) or excluding (2-) the correct quantitation equivalent are also available.

See the next section for additional details about how matching. The search is performed with the `searchGrid` method, possibly on a subset of the data (see Methods and Examples sections for further details).

The parameters used for matching can be set manually with `setPpmError`, `setRtNsd`, `setImDiff` respectively, or using `setBestGridParams` to apply best parameters as defined using the grid search. See example and method documentation for details.

**Identification transfer: matching identification peptides and quantitation EMRTs:** The identification peptide - quantitation EMRT matching, termed identification transfer, is performed using the best parameters, as defined above with a grid search, or using user-defined parameters. Matching is considered successful when one and only one EMRT is found in the mass tolerance/retention time/ion mobility window defined by the error ppm, number of retention time standard deviations, and ion mobility difference parameters. The values of uniquely matched EMRTs are reported in the final *matched* dataframe that can be exported (see below). If however, none or more than one EMRTs are matched, 0 or the number of matches are reported.

As identification peptides are serially individually matched to 'close' EMRTs, it is possible for peptides to be matched the same EMRT independently. Such cases are reported as -1 in the results dataframes.

The results can be assessed using the `plotEMRTtable` (or `getEMRTtable` to retrieve the values) and performance methods. The former shows the number of identification peptides assigned to none (0), exactly 1 (1) or more ( $> 2$ ) EMRTs. The latter method reports matched identification peptides, the number of (q-value and protein FPR filtered) identification and quantitation peptides. Matched EMRT and quantitation peptide numbers are then compared calculating the synapter enrichment ( $100 * (\text{synapter} - \text{quant}) / \text{quant}$ ) and Venn counts.

**Remove Less Intense Peaks:** As an additional step it is possible to remove less intense peaks from the spectra and fragment data. Use `plotCumulativeNumberOfFragments` to plot the number of fragments vs the intensity and to find a good threshold. The `filterFragments` method could remove peaks if the intensity is below a specified threshold via the `minIntensity` argument. Set the `maxNumber` argument to keep only the `maxNumber` highest peaks/fragments. The `what` argument controls the data on which the filter is applied. Use `what = "fragments.ident"` for the identification fragments and `what = "spectra.quant"` for the quantitation spectra data.

**Fragment Matching:** After importing fragment and spectra data it is possible to match peaks between the identification fragments and the quantitation spectra using the `fragmentMatching` method. Use `setFragmentMatchingPpmTolerance` to set the maximal allowed tolerance for considering a peak as identical. There are two different methods to visualise the results of the fragment matching procedure. `plotFragmentMatching` plots the fragments and spectra for each considered pair. `plotFragmentMatchingPerformance` draws two plots. On the left panel you could see the performance of different thresholds for the number of common peaks for unique matches. The right panel visualizes the performance of different differences (delta) of common peaks between the best match (highest number of common peaks) and the second best match in each non unique match group. `plotFragmentMatchingPerformance` returns the corresponding values invisible or use `fragmentMatchingPerformance` to access these data. Use `filterUniqueMatches` and `filterNonUniqueMatches` to remove unique or non unique matches below the threshold of common peaks respective the difference in common peaks from the `MatchedEMRTs` data.frame.

**Exporting and saving data:** The merged identification and quantitation peptides can be exported to csv using the `writeMergedPeptides` method. Similarly, the matched identification peptides and quantitation EMRTs are exported with `writeMatchedEMRTs`.

Complete Synapter instances can be serialised with `save`, as any R object, and reloaded with `load` for further analysis.

It is possible to get the fragment and spectra data from the identification and quantitation run using `getIdentificationFragments` respectively `getQuantitationSpectra`.

## Methods

### Analysis methods:

**mergePeptides** `signature(object = "Synapter")`: Merges quantitation and identification final peptide data, used to perform retention time modelling (see `modelRt` below).

**modelRt** `signature(object = "Synapter", span = "numeric")`: Performs local polynomial regression fitting (see [loess](#)) retention time alignment using span parameter value to control the degree of smoothing.

**modelIntensity** `signature(object = "Synapter", span = "numeric")`: Performs local polynomial regression fitting (see [loess](#)) intensity values using span parameter value to control the degree of smoothing.

**findEMRTs** signature(object = "Synapter", ppm = "numeric", nsd = "numeric", imdiff = "numeric"): Finds EMRTs matching identification peptides using ppm mass tolerance, nsd number of retention time standard deviations and imdiff difference in ion mobility. The last three parameters are optional if previously set with setPpmError, setRtNsd, setImDiff, or, better, setBestGridParams (see below).

**rescueEMRTs** signature(object = "Synapter", method = c("rescue", "copy")): The method parameter defined the behaviour for those high confidence features that were identified in both identification and quantitation acquisitions and used for the retention time model (see mergePeptides). Prior to version 1.1.1, these features were transferred from the quantitation pep3d file if unique matches were found, as any feature ("transfer"). As a result, those matching 0 or > 1 EMRTs were quantified as NA. The default is now to "rescue" the quantitation values of these by directly retrieving the data from the quantification peptide data. Alternatively, the quantitation values for these features can be directly taken from the quantitation peptide data using "copy", thus effectively bypassing identification transfer.

**searchGrid** signature(object="Synapter", ppms="numeric", nsds="numeric", imdiffs = "numeric", subset="numeric", n = "numeric", verbose="logical"): Performs a grid search. The grid is defined by the ppm, nsd and imdiffs numerical vectors, representing the sequence of values to be tested. Default are seq(5, 20, 2), seq(0.5, 5, 0.5), seq(0.2, 2, 0.2) respectively. To ignore ion mobility set imdiffs = Inf. subset and n allow to use a randomly chosen subset of the data for the grid search to reduce search time. subset is a numeric, between 0 and 1, describing the percentage of data to be used; n specifies the absolute number of feature to use. The default is to use all data. verbose controls whether textual output should be printed to the console. (Note, the mergedEMRTs value used in internal calls to findEMRTs is "transfer" - see findEMRTs for details).

**fragmentMatching** signature(object="Synapter", ppm = "numeric", verbose = "logical"): Performs a fragment matching between spectra and fragment data. The ppm argument controls the tolerance that is used to consider two peaks (MZ values) as identical. If verbose is TRUE (default) a progress bar is shown.

*Methods to display, access and set data:*

**show** signature(object = "Synapter"): Display object by printing a summary to the console.

**dim** signature(x="Synapter"): Returns a list of dimensions for the identification peptide, quantitation peptide, merged peptides and matched features data sets.

**inputFiles** signature(object="Synapter"): Returns a character of length 6 with the names of the input files used as identpeptide, quantpeptide, quantpep3d, fasta, identfragments and quantspectra.

**getLog** signature(object="Synapter"): Returns a character of variable length with a summary of processing undergone by object.

**getGrid** signature(object="Synapter", digits = "numeric"): Returns a named list of length 3 with the percent of total (prcntTotal), percent of model (prcntModel) and detailed (details) grid search results. The details grid search reports the proportion of correctly assigned features (+1) to all unique assignments (+1 and -1). Values are rounded to 3 digits by default.

**getGridDetails** signature(object="Synapter"): Returns a list of number of ..., -2, -1, 0, +1, +2, ... results found for each of the ppm/nsd pairs tested during the grid search.

**getBestGridValue** signature(object="Synapter"): Returns a named numeric of length 3 with best grid values for the 3 searches. Names are prcntTotal, prcntModel and details.

**getBestGridParams** signature(object="Synapter"): Returns a named list of matrices (prcntTotal, prcntModel and details). Each matrix gives the ppm and nsd pairs that yielded the best grid values (see getBestGridValue above).

**setBestGridParams** signature(object="Synapter", what="character"): This methods set the best parameter pair, as determined by what. Possible values are auto (default), model, total and details. The 3 last ones use the (first) best parameter values as reported by getBestGridParams. auto uses the best model parameters and, if several best pairs exists, the one that maximises details is selected.

**setPepScoreFdr** signature(object="Synapter", fdr = "numeric"): Sets the peptide score false discovery rate (default is 0.01) threshold used by filterQuantPepScore and filterIdentPepScore.

**getPepScoreFdr** signature(object="Synapter"): Returns the peptide false discovery rate threshold.

**setIdentPpmError** signature(object="Synapter", ppm = "numeric"): Set the identification mass tolerance to ppm (default 10).

**getIdentPpmError** signature(object="Synapter"): Returns the identification mass tolerance.

**setQuantPpmError** signature(object="Synapter", ppm = "numeric"): Set the quantitation mass tolerance to ppm (default 10).

**getQuantPpmError** signature(object="Synapter"): Returns the quantitation mass tolerance.

**setPpmError** signature(object="Synapter", ppm = "numeric"): Sets the identification and quantitation mass tolerance ppm (default is 10).

**setLowessSpan** signature(object="Synapter", span = "numeric"): Sets the loess span parameter; default is 0.05.

**getLowessSpan** signature(object="Synapter"): Returns the span parameter value.

**setRtNsd** signature(object="Synapter", nsd = "numeric"): Sets the retention time tolerance nsd, default is 2.

**getRtNsd** signature(object="Synapter"): Returns the value of the retention time tolerance nsd.

**setImDiff** signature(object="Synapter", imdiff = "numeric"): Sets the ion mobility tolerance imdiff, default is 0.5.

**getImDiff** signature(object="Synapter"): Returns the value of the ion mobility tolerance imdiff.

**getPpmErrorQs** signature(object="Synapter", qs = "numeric", digits = "numeric"): Returns the mass tolerance qs quantiles (default is c(0.25, 0.5, 0.75, seq(0.9, 1, 0.01)) for the identification and quantitation peptides. Default is 3 digits.

**getRtQs** signature(object="Synapter", qs = "numeric", digits = "numeric"): Returns the retention time tolerance qs quantiles (default is c(0.25, 0.5, 0.75, seq(0.9, 1, 0.01)) for the identification and quantitation peptides. Default is 3 digits.

**getPepNumbers** signature(object="Synapter"): Returns the number of regular and random quantitation and identification peptide considered for p-value calculation and used to plot the score densities (see plotPepScores). Especially the difference between random and regular entries are informative in respect with the confidence of the random scores distribution.

**setFragmentMatchingPpmTolerance** signature(object="Synapter", ppm = "numeric"): Sets the fragment matching mass tolerance ppm (default is 25).

**getFragmentMatchingPpmTolerance** signature(object="Synapter"): Returns the fragment matching mass tolerance in ppm.

**showFdrStats** signature(object="Synapter", k = "numeric"): Returns a named list of length 2 with the proportion of identification and quantitation peptides that are considered significant with a threshold of k (default is c(0.001, 0.01, 0.5, 0.1)) using raw and adjusted p-values/q-values.

**getEMRTtable** signature(object="Synapter"): Returns a table with the number of 0, 1, 2, ... assigned EMRTs.

**performance** signature(object="Synapter", verbose = TRUE): Returns (and displays, if verbose) the performance of the synapter analysis.

**performance2** signature(object="Synapter", verbose = TRUE): Returns (and displays, if verbose) information about number of missing values and identification source of transferred EMRTs.

**fragmentMatchingPerformance** signature(object="Synapter", what = c("unique", "non-unique")): Returns the performance of the fragment matching for unique or non-unique matches. The return value is a matrix with seven columns. The first column ncommon/deltacommon contains the thresholds. Column 2 to 5 are the true positives tp, false positives fp, true negatives tn, false negatives fn for the merged peptide data. The sixth column all shows the corresponding number of peptides for all peptides (not just the merged ones) and the last column shows the FDR fdr for the current threshold (in that row) for the merged data. Please note that the FDR is overfitted/underestimated because the merged peptides are the peptides from the highest quality spectra where PLGS could easily identify the peptides. The peptides that are not present in the merged data are often of lower quality hence the FDR would be higher by trend.  
See plotFragmentMatchingPerformance for a graphical representation.

#### Filters:

**filterUniqueDbPeptides** signature(object="Synapter", missedCleavages = 0, IisL = TRUE, verbose = TRUE): This method first digests the fasta database file and keeps unique tryptic peptides. (NOTE: since version 1.5.3, the tryptic digestion uses the cleaver package, replacing the more simplistic inbuild function. The effect of this change is documented in <https://github.com/lgatto/synapter/pull/47>). The number of maximal missed cleavages can be set as missedCleavages (default is 0). If IisL = TRUE Isoleucin and Leucin are treated as the same aminoacid. In this case sequences like "ABCI", "ABCL" are removed because they are not unique anymore. If IisL = FALSE (default) "ABCI" and "ABCL" are reported as unique. The peptide sequences are then used as a filter against the identification and quantitation peptides, where only unique proteotypic instances (no miscleavage allowed by default) are eventually kept in the object instance. This method also removes any additional duplicated peptides, that would not match any peptides identified in the fasta database.

**filterUniqueQuantDbPeptides** signature(object="Synapter", missedCleavages = 0, IisL = TRUE, verbose = TRUE): As filterUniqueDbPeptides for quantitation peptides only.

**filterUniqueIdentDbPeptides** signature(object="Synapter", missedCleavages = 0, IisL = TRUE, verbose = TRUE): As filterUniqueDbPeptides for identification peptides only.

**filterQuantPepScore** signature(object="Synapter", fdr = "numeric", method = "character"): Filters the quantitation peptides using fdr false discovery rate. fdr is missing by default and is retrieved with getPepScoreFdr automatically. If not set, default value of 0.01 is

used. method defines how to perform p-value adjustment; one of BH, Bonferrone or qval.  
See details section for more information.

**filterIdentPepScore** signature(object="Synapter", fdr = "numeric", method = "character"): As filterQuantPepScore, but for identification peptides.

**filterQuantProtFpr** signature(object="Synapter", fpr = "numeric"): Filters quantitation peptides using the protein false positive rate (erroneously defined as a FPR, should be FDR), as reported by PLGS, using threshold set by fpr (missing by default) or retrieved by getProtFpr.

**filterIdentProtFpr** signature(object="Synapter", fpr = "numeric"): as filterQuantProtFpr, but for identification peptides.

**filterQuantPpmError** signature(object="Synapter", ppm = "numeric"): Filters the quantitation peptides based on the mass tolerance ppm (default missing) provided or retrieved automatically using getPpmError.

**filterIdentPpmError** signature(object="Synapter"): as filterQuantPpmError, but for identification peptides.

**filterFragments** signature(object = "Synapter", what = c("fragments.ident", "spectra.quant"), minIntensity = "numeric", maxNumber = "numeric", verbose = "logical"): Filters the spectra/fragment data using a minimal intensity threshold (minIntensity) or a maximal number of peaks/fragments threshold (maxNumber). Please note that the maximal number is transferred to an intensity threshold and the result could contain less peaks than specified by maxNumber. If both arguments are given, the more aggressive one is chosen. Use the what argument to specify the data that should be filtered. Set what = "fragments.ident" for the identification fragment data or what = "spectra.quant" for the quantitation spectra. If verbose is TRUE (default) a progress bar is shown.

**filterUniqueMatches** signature(object="Synapter", minNumber = "numeric"): Removes all unique matches that have less than minNumber of peaks/fragments in common. Use fragmentMatchingPerformance(..., what="unique")/plotFragmentMatchingPerformance (left panel) to find an ideal threshold.

**filterNonUniqueMatches** signature(object="Synapter", minDelta = "numeric"): Removes all non unique matches that have a difference between the best match (highest number of common peaks/fragments, treated as true match) and the second best match (second highest number of common peaks/fragments) less than minDelta. For the matches above the threshold only the one with the highest number of common peaks/fragments in each match group is kept. Use fragmentMatchingPerformance(..., what="non-unique")/plotFragmentMatchingPerformance (right panel) to find an ideal threshold.

**filterNonUniqueIdentMatches** signature(object="Synapter"): Removes all non unique identification matches. In rare circumstances (if the grid search parameters are too wide/relaxed or a fragment library is used) it could happen that the searchGrid methods matches a single quantification EMRT to multiple identification EMRTs. This method removes all these non unique matches.

*Plotting:*

**plotPpmError** signature(object="Synapter", what = "character"): Plots the proportion of data against the mass error tolerance in ppms. Depending on what, the data for identification (what = "Ident"), quantitation (what = "Quant") or "both" is plotted.

**plotRtDiffs** signature(object="Synapter", ...): Plots a histogram of retention time differences after alignments. ... is passed to hist.

**plotRt** signature(object="Synapter", what = "character", f = "numeric", nsd = "numeric"):

Plots the Identification - Quantitation retention time difference as a function of the Identification retention time. If what is "data", two plots are generated: one ranging the full range of retention time differences and one focusing on the highest data point density and showing models with various span parameter values, as defined by f (default is 2/3, 1/2, 1/4, 1/10, 1/16, 1/25, 1/50, passed as a numeric). If what is "model", a focused plot with the applied span parameter is plotted and areas of nsd (default is x(1, 3, 5) number of standard deviations are shaded around the model.

**plotIntensity** signature(object="Synapter", what = "character", f = "numeric", nsd = "numeric"):

Plots the (log2) ratio of Identification and Quantitation intensities as a function of the Identification retention time. If what is "data", two plots are generated: one ranging the full range of ratios and one focusing on the highest data point density and showing models with various span parameter values, as defined by f (default is 2/3, 1/2, 1/4, 1/10, 1/16, 1/25, 1/50, passed as a numeric). If what is "model", a focused plot with the applied span parameter is plotted and areas of nsd (default is x(1, 3, 5) number of standard deviations are shaded around the model.

**plotPepScores** signature(object="Synapter"): Plots the distribution of random and regular peptide scores for identification and quantitation features. This reflects how peptide p-values are computed. See also getPepNumbers.

**plotFdr** signature(object="Synapter", method = "character"): Displays 2 plots per identification and quantitation peptides, showing the number of significant peptides as a function of the FDR cut-off and the expected false number of false positive as a number of significant tests. PepFrag 1 and 2 peptides are illustrated on the same figures. These figures are adapted from [plot.qvalue](#). method, one of "BH", "Bonferroni" or "qval", defines what identification statistics to use.

**plotEMRTtable** signature(object="Synapter"): Plots the barchart of number of 0, 1, 2, ... assigned EMRTs (see getEMRTtable).

**plotGrid** signature(object="Synapter", what = "character"), maindim = "character"):

Plots a heatmap of the respective grid search results. This grid to be plotted is controlled by what: "total", "model" or "details" are available. If ion mobility was used in the grid search you can use maindim to decide which dimensions should be shown. maindim could be one of "im" (default), "rt" and "mz". If maindim = "im" a heatmap for each ion mobility threshold is drawn. For maindim = "rt" and maindim you get a heatmap for each retention time respective mass threshold.

**plotFeatures** signature(object="Synapter", what = "character", xlim = "numeric", ylim = "numeric", ionmobility = "logical"):

Plots the retention time against precursor mass space. If what is "all", three (six if ion mobility is available and ionmobility = TRUE (default is FALSE); three additional plots with precursor mass against ion mobility) such plots are created side by side: for the identification peptides, the quantitation peptides and the quantitation Pep3D data. If what is "some", a subset of the rt/mass space can be defined with xlim (default is c(40, 60)) and ylim (default is c(1160, 1165)) and identification peptide, quantitation peptides and EMRTs are presented on the same graph as grey dots, blue dots and red crosses respectively. In addition, rectangles based on the ppm and nsd defined tolerances (see setPpmError and setNsdError) are drawn and centered at the expected modelled retention time. This last figure allows to visualise the EMRT matching.

**plotFragmentMatching** signature(object = "Synapter", key = "character", column = "character",

verbose = "logical", ...): Plots two spectra and fragments against each other. Please see [plotFragmentMatching](#) for details.

**plotFragmentMatchingPerformance** signature(object = "Synapter", showAllPeptides = FALSE): Creates two plots. The left panel shows the performance of filtering the unique matches of the merged peptides using a different number of common peaks. The right panel shows the performance of filtering the non unique matches of the merged peptides using different differences (delta) in common peaks/fragments. These differences (delta) are calculated between the match with the highest number of common peaks/fragments and the second highest one. Use `filterUniqueMatches` and `filterNonUniqueMatches` to filter the MatchedEMRT data.frame using one of these thresholds. This function returns a list with two named elements (unique and nonunique invisibly. These are the same data as return by `fragmentMatchingPerformance`.

Use `showAllPeptides=TRUE` to add a line for all peptides (not just the merged ones) to both plots.

**plotCumulativeNumberOfFragments** signature(object = "Synapter", what = c("fragments.ident", "spectra.quant")): Plots the cumulative number of the fragments/peaks vs their intensity (log10 scaled). Use the `what` argument to create this plot for the identification fragments (what = "fragments.quant") or the the quantitation spectra (what = "spectra.quant").

*Exporters:*

**writeMergedPeptides** signature(object="Synapter", file = "character", what = "character", ...): Exports the merged peptide data to a comma-separated file (default name is "Res-MergedPeptides.csv").

**writeMatchedEMRTs** signature(object="Synapter", file = "character", ...): As above, saving the matched EMRT table.

**writeIdentPeptides** signature(object="Synapter", file = "character", ...): As above, exporting the identification peptide data.

**writeQuantPeptides** signature(object="Synapter", file = "character", ...): A above, exporting the quantitation peptide data.

**getIdentificationFragments** signature(object="Synapter"): returns the identification fragments as [MSnExp](#).

**getQuantitationSpectra** signature(object="Synapter"): returns the quantitation spectra as [MSnExp](#).

*Other:*

**as(, "MSnSet")** signature(x = "Synapter"): Coerce object from Synapter to MSnSet class.

**validObject** signature(object = "Synapter"): Test whether a given Synapter object is valid.

**updateObject** signature(object = "Synapter"): Updates an old Synapter object.

## Author(s)

Laurent Gatto <lg390@cam.ac.uk>

## References

Käll L, Storey JD, MacCoss MJ, Noble WS Posterior error probabilities and false discovery rates: two sides of the same coin. *J Proteome Res.* 2008a Jan; 7:(1)40-4

Bonferroni single-step adjusted p-values for strong control of the FWER.

Benjamini Y. and Hochberg Y. Controlling the false discovery rate: a practical and powerful approach to multiple testing. *J. R. Statist. Soc. B.*, 1995, Vol. 57: 289-300.

Storey JD and Tibshirani R. Statistical significance for genome-wide experiments. *Proceedings of the National Academy of Sciences*, 2003, 100: 9440-9445.

Käll, Storey JD, MacCoss MJ, Noble WS Assigning significance to peptides identified by tandem mass spectrometry using decoy databases. *J Proteome Res*. 2008b Jan; 7:(1)29-34

Improving qualitative and quantitative performance for MSE-based label free proteomics, N.J. Bond, P.V. Shliaha, K.S. Lilley and L. Gatto, *Journal of Proteome Research*, 2013, in press.

The Effects of Travelling Wave Ion Mobility Separation on Data Independent Acquisition in Proteomics Studies, P.V. Shliaha, N.J. Bond, L. Gatto and K.S. Lilley, *Journal of Proteome Research*, 2013, in press.

Trypsin cleavage:

Glatter, Timo, et al. Large-scale quantitative assessment of different in-solution protein digestion protocols reveals superior cleavage efficiency of tandem Lys-C/trypsin proteolysis over trypsin digestion. *Journal of proteome research* 11.11 (2012): 5145-5156. <http://dx.doi.org/10.1021/pr300273g>

Rodriguez, Jesse, et al. Does trypsin cut before proline?. *Journal of proteome research* 7.01 (2007): 300-305. <http://dx.doi.org/10.1021/pr0705035>

Brownridge, Philip, and Robert J. Beynon. The importance of the digest: proteolysis and absolute quantification in proteomics. *Methods* 54.4 (2011): 351-360. <http://dx.doi.org/10.1016/j.jymeth.2011.05.005>

cleaver's rules are taken from: [http://web.expasy.org/peptide\\_cutter/peptidecutter\\_enzymes.html#Tryps](http://web.expasy.org/peptide_cutter/peptidecutter_enzymes.html#Tryps)

## Examples

```
library(synapter) ## always needed

## Not run:
## (1) Construction - to create your own data objects
synapterTiny <- Synapter()

## End(Not run)

## let's use synapterTiny, shipped with the package
synapterTinyData() ## loads/prepares the data
synapterTiny ## show object

## (2) Filtering
## (2.1) Peptide scores and FDR

## visualise/explore peptide id scores
plotPepScores(synapterTiny)
getPepNumbers(synapterTiny)

## filter data
filterUniqueDbPeptides(synapterTiny) ## keeps unique proteotypic peptides
filterPeptideLength(synapterTiny, l = 7) ## default length is 7

## visualise before FDR filtering
```

```
plotFdr(synapterTiny)

setPepScoreFdr(synapterTiny, fdr = 0.01) ## optional
filterQuantPepScore(synapterTiny, fdr = 0.01) ## specifying FDR
filterIdentPepScore(synapterTiny) ## FDR not specified, using previously set value

## (2.2) Mass tolerance
getPpmErrorQs(synapterTiny)
plotPpmError(synapterTiny, what="Ident")
plotPpmError(synapterTiny, what="Quant")

setIdentPpmError(synapterTiny, ppm = 20) ## optional
filterQuantPpmError(synapterTiny, ppm = 20)
## setQuantPpmError(synapterTiny, ppm = 20) ## set quant ppm threshold below
filterIdentPpmError(synapterTiny, ppm=20)

filterIdentProtFpr(synapterTiny, fpr = 0.01)
filterQuantProtFpr(synapterTiny, fpr = 0.01)

getPpmErrorQs(synapterTiny) ## to be compared with previous output

## (3) Merge peptide sequences
mergePeptides(synapterTiny)

## (4) Retention time modelling
plotRt(synapterTiny, what="data")
setLowessSpan(synapterTiny, 0.05)
modelRt(synapterTiny) ## the actual modelling
getRtQs(synapterTiny)
plotRtDiffs(synapterTiny)
## plotRtDiffs(synapterTiny, xlim=c(-1, 1), breaks=500) ## pass parameters to hist()
plotRt(synapterTiny, what="model") ## using default nsd 1, 3, 5
plotRt(synapterTiny, what="model", nsd=0.5) ## better focus on model

plotFeatures(synapterTiny, what="all")
setRtNsd(synapterTiny, 3) ## RtNsd and PpmError are used for detailed plot
setPpmError(synapterTiny, 10) ## if not set manually, default values are set automatically
plotFeatures(synapterTiny, what="some", xlim=c(36,44), ylim=c(1161.4, 1161.7))
## best plotting to svg for zooming

set.seed(1) ## only for reproducibility of this example

## (5) Grid search to optimise EMRT matching parameters
searchGrid(synapterTiny,
           ppms = 7:10, ## default values are 5, 7, ..., 20
           nsds = 1:3, ## default values are 0.5, 1, ..., 5
           subset = 0.2) ## default is 1
## alternatively, use 'n = 1000' to use exactly
## 1000 randomly selected features for the grid search
getGrid(synapterTiny) ## print the grid
getGridDetails(synapterTiny) ## grid details
plotGrid(synapterTiny, what = "total") ## plot the grid for total matching
plotGrid(synapterTiny, what = "model") ## plot the grid for matched modelled feature
```

```

plotGrid(synapterTiny, what = "details") ## plot the detail grid
getBestGridValue(synapterTiny) ## return best grid values
getBestGridParams(synapterTiny) ## return parameters corresponding to best values
setBestGridParams(synapterTiny, what = "auto") ## sets RtNsd and PpmError according the grid results
## 'what' could also be "model", "total" or "details"
## setPpmError(synapterTiny, 12) ## to manually set values
## setRtNsd(synapterTiny, 2.5)

## (6) Matching ident peptides and quant EMRTs
findEMRTs(synapterTiny)
plotEMRTtable(synapterTiny)
getEMRTtable(synapterTiny)
performance(synapterTiny)
performance2(synapterTiny)

## (7) Exporting data to csv spreadsheets
writeMergedPeptides(synapterTiny)
writeMergedPeptides(synapterTiny, file = "myresults.csv")
writeMatchedEMRTs(synapterTiny)
writeMatchedEMRTs(synapterTiny, file = "myresults2.csv")
## These will export the filter peptide data
writeIdentPeptides(synapterTiny, file = "myIdentPeptides.csv")
writeQuantPeptides(synapterTiny, file = "myQuantPeptides.csv")
## If used right after loading, the non-filted data will be exported

```

---

synapterGuide

*Opens the 'synapter' vignette*


---

## Description

Opens the relevant vignette; a shortcut to using vignette. `synapterGuide()` gives access to the main overview vignette.

## Usage

```
synapterGuide()
```

## Author(s)

Laurent Gatto

---

```
synapterPlgsAgreement, MSnSet-method
```

*Synapter/PLGS Agreement*

---

## Description

This method checks the agreement between synapter analysis and PLGS results.

## Usage

```
## S4 method for signature 'MSnSet'
synapterPlgsAgreement(object, ...)
```

## Arguments

object	An <a href="#">MSnSet</a> object.
...	further arguments, not used yet.

## Details

Each synapter object has synapterPlgsAgreement column in its MatchedEMRTs data.frame (see [writeMatchedEMRTs](#)). After converting the synapter object into an [MSnSet](#) instance via `as(synapterobject, "MSnSet")` this column could be find in the feature data (`fData(msnset)$synapterPlgsAgreement`).

In the synapterPlgsAgreement each peptide is classified as:

- "agree": EMRT identified in identification and quantitation run by PLGS and same EMRT matched in synapter's grid search.
- "disagree": EMRT identified in identification and quantitation run by PLGS and a different EMRT was matched in synapter's grid search.
- "no\_plgs\_id": EMRT was *not* identified in the quantitation run by PLGS but matched in synapter's grid search.
- "no\_synapter\_transfer": EMRT was identified in the identification and quantitation run by PLGS but not matched in synapter's grid search.
- "no\_id\_or\_transfer": EMRT was *not* identified in the quantitation run by PLGS and not matched in synapter's grid search.
- "multiple\_ident\_matches": a single quantitation EMRT was matched by synapter to multiple identification EMRTs found by PLGS (could happen if the grid search parameters are too relaxed).

After combining multiple [MSnSet](#) the method synapterPlgsAgreement adds additional columns to the feature data:

- nIdentified: how often a peptide was identified across multiple runs?
- nAgree: how often a peptide was identified by PLGS and synapter across multiple runs (counts "agree" entries)?
- nDisagree: how often a peptide was differently identified by PLGS and synapter across multiple runs (counts "disagree" entries)?
- synapterPlgsAgreementRatio:  $nAgree / (nAgree + nDisagree)$ .

**Value**

[MSnSet](#) where the columns nIdentified, nAgree, nDisagree and synapterPlgsAgreementRatio were added to the feature data.

**Author(s)**

Sebastian Gibb <mail@sebastiangibb.de>

**References**

See discussion on github: <https://github.com/lgatto/synapter/issues/73>

**See Also**

MSnSet documentation: [MSnSet](#)

---

synapterTinyData	<i>Loads a small test data for the 'synapter' package</i>
------------------	---

---

**Description**

Instead of using data to load the synapterTiny data set, synapterTinyData will load it and initialise it for proper downstream analysis, during which the 04\_test\_database.fasta file, provided with the package and references inside the object needs, to be accessed. However, as the exact location can not be known in advance, the reference is updated with the file's correct local path.

This data set has been generated with the [Synapter](#) constructor. Note that the input data file sizes have been reduced by depleting many rows (peptides and EMRTs) from the original csv files.

In addition, several columns that were not necessary for processing were also removed. As such, the data stored in synapterTiny does not reflect the data obtained when following the section 'Preparing the input data' in the section, without however affecting the processing and final results.

**Usage**

```
synapterTinyData()
```

**Value**

A character vector with the data set name, "synapterTiny". Used for its side effect of loading synapterTiny, an instance of class [Synapter](#), in .GlobalEnv.

**Author(s)**

Laurent Gatto

## Source

Improving qualitative and quantitative performance for MSE-based label free proteomics, N.J. Bond, P.V. Shliaha, K.S. Lilley and L. Gatto, Journal of Proteome Research, 2013, in press.

The Effects of Travelling Wave Ion Mobility Separation on Data Independent Acquisition in Proteomics Studies, P.V. Shliaha, N.J. Bond, L. Gatto and K.S. Lilley, Journal of Proteome Research, 2013, in press.

## Examples

```
synapterTinyData()
synapterTiny
```

---

synergise	<i>Synergise identification and quantitation results</i>
-----------	--

---

## Description

Performs a complete default analysis on the files defined in filenames, creates a complete html report and saves/exports all results as csv and rds files. See details for a description of the pipeline and [Synapter](#) for manual execution of individual steps.

## Usage

```
synergise(filenames, master = FALSE, object, outputdir,
  outputfile = paste0("synapter_report_", strftime(Sys.time(),
    "%Y%m%d-%H%M%S"), ".html"), fdr = 0.01, fdrMethod = c("BH",
    "Bonferroni", "qval"), fpr = 0.01, peplen = 7, missedCleavages = 0,
    IisL = FALSE, identppm = 20, quantppm = 20, uniquepep = TRUE,
    span = 0.05, grid.ppm.from = 2, grid.ppm.to = 20,
    grid.ppm.by = 2, grid.nsd.from = 0.5, grid.nsd.to = 5,
    grid.nsd.by = 0.5, grid.subset = 1, grid.n = 0,
    grid.param.sel = c("auto", "model", "total", "details"),
    mergedEMRTs = c("rescue", "copy", "transfer"),
    template = system.file("reports", "synergise1.Rmd", package =
    "synapter"), verbose = interactive())
```

## Arguments

filenames	A named list of file names to be load. The names must be identpeptide, quantpeptide, quantpep3d and fasta (could be an RDS file created by <code>link{createUniquePeptideDb}</code> identpeptide can be a csv final peptide file (from PLGS) or a saved " <a href="#">MasterPeptides</a> " data object as created by <a href="#">makeMaster</a> if working with <i>master</i> peptide data. To serialise the " <a href="#">MasterPeptides</a> " instance, use the <code>saveRDS</code> function, and file extension rds.
master	A logical indicating if the identification final peptide files are master (see <a href="#">makeMaster</a> ) or regular files. Default is FALSE.

object	An instance of class <i>Synapter</i> that will be copied, processed and returned. If filenames are also provided, the latter and object's <code>inputFiles</code> will be checked for equality.
outputdir	A character with the full path to an existing directory.
outputfile	A character with the file name for the report.
fdr	Peptide false discovery rate. Default is 0.01.
fdrMethod	P-value adjustment method. One of "BH" (default) for Benjamini and Hochberg (1995), "Bonferroni" for Bonferroni's single-step adjusted p-values for strong control of the FWER and "qval" from the <i>qvalue</i> package. See <a href="#">Synapter</a> for references.
fpr	Protein false positive rate. Default is 0.01.
peplen	Minimum peptide length. Default is 7.
missedCleavages	Number of allowed missed cleavages. Default is 0.
IisL	If TRUE Isoleucin and Leucin are treated as equal. In this case sequences like "ABCI", "ABCL" are removed because they are not unique. If FALSE (default) "ABCI" and "ABCL" are reported as unique.
identppm	Identification mass tolerance (in ppm). Default is 20.
quantppm	Quantitation mass tolerance (in ppm). Default is 20.
uniquepep	A logical is length 1 indicating if only unique peptides in the identification and quantitation peptides as well as unique tryptic peptides as defined in the fasta file. Default is TRUE.
span	The loess span parameter. Default is 0.05.
grid.ppm.from	Mass tolerance (ppm) grid starting value. Default is 2.
grid.ppm.to	Mass tolerance (ppm) grid ending value. Default is 20.
grid.ppm.by	Mass tolerance (ppm) grid step value. Default is 2.
grid.nsd.from	Number of retention time stdev grid starting value. Default is 0.5.
grid.nsd.to	Number of retention time stdev grid ending value. Default is 5.
grid.nsd.by	Number of retention time stdev grid step value. Default is 0.5.
grid.subset	Percentage of features to be used for the grid search. Default is 1.
grid.n	Absolute number of features to be used for the grid search. Default is 0, i.e ignored.
grid.param.sel	Grid parameter selection method. One of auto (default), details, model or total. See <a href="#">Synapter</a> for details on these selection methods.
mergedEMRTs	One of "rescue" (default), "copy" or "transfer". See the documentation for the <code>findEMRTs</code> function in <a href="#">Synapter</a> for details.
template	A character full path to Rmd template.
verbose	A logical indicating if progress output should be printed to the console. Default is TRUE.

## Details

Data can be input as a [Synapter](#) object if available or as a list of files (see `filenames`) that will be used to read the data in. The html report and result files will be created in the `outputdir` folder. All other input parameters have default values.

The data processing and analysis pipeline is as follows:

1. If `uniquepep` is set to `TRUE` (default), only unique proteotypic identification and quantitation peptides are retained.
2. Peptides are filtered for a  $FDR \leq fdr$  (default is 0.01) using the "BH" method (see `fdr` and `fdrMethod` parameters for details).
3. Peptide with a mass tolerance  $> 20$  ppm (see `quantppm` and `identppm`) are filtered out.
4. Peptides with a protein false positive rate (as reported by the PLGS software)  $> fpr$  are filtered out.
5. Common identification and quantitation peptides are merged and a retention time model is created using the Local Polynomial Regression Fitting ([loess](#) function for the `stats` package) using a default span value of 0.05.
6. A grid search to optimise the width in retention time and mass tolerance for EMRTs matching is performed. The default grid search space is from 0.5 to 5 by 0.5 retention time model standard deviations (see `grid.nsd.from`, `grid.nsd.to` and `grid.nsd.by` parameters) and from 2 to 20 by 2 parts per million (ppm) for mass tolerance (see `grid.ppm.from`, `grid.ppm.to` and `grid.ppm.by` parameters). The data can be subset using using an absolute number of features (see `grid.n`) or a fixed percentage (see `grid.subset`). The pair of optimal `nsd` and `ppm` is chosen (see `grid.param.sel` parameter).
7. The quantitation EMRTs are matched using the optimised parameters.

If a master identification file is used (`master` is set to `TRUE`, default is `FALSE`), the relevant actions that have already been executed when the file was created with [makeMaster](#) are not repeated here.

## Value

Invisibly returns an object of class `Synapter`. Used for its side effect of creating an html report of the run in `outputdir`.

## Author(s)

Laurent Gatto, Sebastian Gibb

## References

Bond N. J., Shliaha P.V., Lilley K.S. and Gatto L. (2013) J. Prot. Research.

## Examples

```
output <- tempdir() ## a temporary directory
synapterTinyData()
synergise(object = synapterTiny, outputdir = output, outputfile = "synapter.html",
          grid.subset = 0.2)
htmlReport <- paste0("file:/// ", file.path(output, "synapter.html")) ## the result report
```

```
## Not run:
browseURL(htmlReport) ## open the report with default browser

## End(Not run)
```

synergise2

*Synergise identification and quantitation results*

## Description

Performs a complete default analysis on the files defined in `filenames`, creates a complete html report and saves/exports all results as csv and rds files. See details for a description of the pipeline and [Synapter](#) for manual execution of individual steps.

## Usage

```
synergise2(filenames, master = FALSE, object, outputdir,
  outputfile = paste0("synapter_report_", strftime(Sys.time(),
    "%Y%m%d-%H%M%S"), ".html"), fdr = 0.01, fdrMethod = c("BH",
    "Bonferroni", "qval"), fpr = 0.01, peplen = 7, missedCleavages = 2,
  lisl = FALSE, identppm = 20, quantppm = 20, uniquepep = TRUE,
  span.rt = 0.05, span.int = 0.05, grid.ppm.from = 2,
  grid.ppm.to = 20, grid.ppm.by = 2, grid.nsd.from = 0.5,
  grid.nsd.to = 5, grid.nsd.by = 0.5, grid.imdiffs.from = 0.6,
  grid.imdiffs.to = 1.6, grid.imdiffs.by = 0.2, grid.subset = 1,
  grid.n = 0, grid.param.sel = c("auto", "model", "total", "details"),
  fm.ppm = 25, fm.ident.minIntensity = 70,
  fm.quant.minIntensity = 70, fm.minCommon = 1, fm.minDelta = 1,
  fm.fdr.unique = 0.05, fm.fdr.nonunique = 0.05,
  mergedEMRTs = c("rescue", "copy", "transfer"),
  template = system.file("reports", "synergise2.Rmd", package =
    "synapter"), verbose = interactive())
```

## Arguments

<code>filenames</code>	A named list of file names to be load. The names must be <code>identpeptide</code> , <code>quantpeptide</code> , <code>quantpep3d</code> and <code>fasta</code> (could be an RDS file created by <code>link{createUniquePeptideDb}</code> ). If <code>fragmentmatching</code> should be used <code>identfragments</code> (could be skipped if a master RDS files is used for <code>identpeptide</code> ) and <code>quantspectra</code> have to be given as well. <code>identpeptide</code> can be a csv final peptide file (from PLGS) or a saved " <a href="#">MasterPeptides</a> " data object as created by <a href="#">makeMaster</a> if working with <i>master</i> peptide data. To serialise the " <a href="#">MasterPeptides</a> " instance, use the <code>saveRDS</code> function, and file extension <code>rds</code> .
<code>master</code>	A logical indicating if the identification final peptide files are master (see <a href="#">makeMaster</a> ) or regular files. Default is <code>FALSE</code> .
<code>object</code>	An instance of class <code>Synapter</code> that will be copied, processed and returned. If <code>filenames</code> are also provided, the latter and <code>object</code> 's <code>inputFiles</code> will be checked for equality.

outputdir	A character with the full path to an existing directory.
outputfile	A character with the file name for the report.
fdr	Peptide false discovery rate. Default is 0.01.
fdrMethod	P-value adjustment method. One of "BH" (default) for Benjamini and HochBerg (1995), "Bonferroni" for Bonferroni's single-step adjusted p-values for strong control of the FWER and "qval" from the qvalue package. See <a href="#">Synapter</a> for references.
fpr	Protein false positive rate. Default is 0.01.
peplen	Minimum peptide length. Default is 7.
missedCleavages	Number of allowed missed cleavages. Default is 2.
IisL	If TRUE Isoleucin and Leucin are treated as equal. In this case sequences like "ABCI", "ABCL" are removed because they are not unique. If FALSE (default) "ABCI" and "ABCL" are reported as unique.
identppm	Identification mass tolerance (in ppm). Default is 20.
quantppm	Quantitation mass tolerance (in ppm). Default is 20.
uniquepep	A logical is length 1 indicating if only unique peptides in the identification and quantitation peptides as well as unique tryptic peptides as defined in the fasta file. Default is TRUE.
span.rt	The loess span parameter for retention time correction. Default is 0.05.
span.int	The loess span parameter for intensity correction. Default is 0.05.
grid.ppm.from	Mass tolerance (ppm) grid starting value. Default is 2.
grid.ppm.to	Mass tolerance (ppm) grid ending value. Default is 20.
grid.ppm.by	Mass tolerance (ppm) grid step value. Default is 2.
grid.nsd.from	Number of retention time stdev grid starting value. Default is 0.5.
grid.nsd.to	Number of retention time stdev grid ending value. Default is 5.
grid.nsd.by	Number of retention time stdev grid step value. Default is 0.5.
grid.imdiffs.from	Ion mobility difference grid starting value. value. Default is 0.6.
grid.imdiffs.to	Ion mobility difference grid ending value. Default is 1.6.
grid.imdiffs.by	Ion mobility difference grid step value. Default is 0.2.
grid.subset	Percentage of features to be used for the grid search. Default is 1.
grid.n	Absolute number of features to be used for the grid search. Default is 0, i.e ignored.
grid.param.sel	Grid parameter selection method. One of auto (default), details, model or total. See <a href="#">Synapter</a> for details on these selection methods.
fm.ppm	Fragment Matching tolerance: Peaks in a range of fm.ppm are considered as identical. Default is 25.

<code>fm.ident.minIntensity</code>	Minimal intensity of a Identification fragment to be not filtered prior to Fragment Matching.
<code>fm.quant.minIntensity</code>	Minimal intensity of a peaks in a Quantitation spectra to be not filtered prior to Fragment Matching.
<code>fm.minCommon</code>	Minimal number of peaks that unique matches need to have in common. Default 1.
<code>fm.minDelta</code>	Minimal difference in number of peaks that non-unique matches need to have to be considered as true match. Default 1.
<code>fm.fdr.unique</code>	Minimal FDR to select <code>fm.minCommon</code> automatically (if both values are given the more restrictive one (that filters more fragments) is used). Default 0.05.
<code>fm.fdr.nonunique</code>	Minimal FDR to select <code>fm.minDelta</code> automatically (if both values are given the more restrictive one (that filters more fragments) is used). Default 0.05.
<code>mergedEMRTs</code>	One of "rescue" (default), "copy" or "transfer". See the documentation for the <code>findEMRTs</code> function in <a href="#">Synapter</a> for details.
<code>template</code>	A character full path to Rmd template.
<code>verbose</code>	A logical indicating if progress output should be printed to the console. Default is TRUE.

## Details

In contrast to [synergise1](#) synergise2 extends the default analysis and offers the following unique features:

- Performing 3D grid search (M/Z, Retention Time, Ion Mobility) for HDMSE data.
- Applying intensity correction.
- Filtering results by fragment matching.

Data can be input as a [Synapter](#) object if available or as a list of files (see `filenames`) that will be used to read the data in. The html report and result files will be created in the `outputdir` folder. All other input parameters have default values.

The data processing and analysis pipeline is as follows:

1. If `uniquepep` is set to TRUE (default), only unique proteotypic identification and quantitation peptides are retained.
2. Peptides are filtered for a  $FDR \leq fdr$  (default is 0.01) using the "BH" method (see `fdr` and `fdrMethod` parameters for details).
3. Peptide with a mass tolerance  $> 20$  ppm (see `quantppm` and `identppm`) are filtered out.
4. Peptides with a protein false positive rate (as reported by the PLGS software)  $> fpr$  are filtered out.
5. Common identification and quantitation peptides are merged and a retention time model is created using the Local Polynomial Regression Fitting ([loess](#) function for the `stats` package) using a default `span.rt` value of 0.05.

6. A grid search to optimise the width in retention time and mass tolerance (and ion mobility for HDMSE) for EMRTs matching is performed. The default grid search space is from 0.5 to 5 by 0.5 retention time model standard deviations (see `grid.nsd.from`, `grid.nsd.to` and `grid.nsd.by` parameters) and from 2 to 20 by 2 parts per million (ppm) for mass tolerance (see `grid.ppm.from`, `grid.ppm.to` and `grid.ppm.by` parameters). If HDMSE data are used the search space is extended from ion mobility difference 0.6 to 1.6 by 0.2 (see `grid.imdiffs.from`, `grid.imdiffs.to` and `grid.imdiffs.by`). The data can be subset using using an absolute number of features (see `grid.n`) or a fixed percentage (see `grid.subset`). The pair of optimal nsd and ppm is chosen (see `grid.param.sel` parameter).
7. Fragment Matching is used to filter false-positive matches from the grid search using a default of 1 common peak for unique matches and at least a difference of 1 common peaks to choose the correct match out of non-unique matches (see `fm.minCommon` and `fm.minDelta`).
8. The intensity is corrected by a Local Polynomial Regression Fitting ([loess](#) function for the `stats` package) using a default `span.int` value of 0.05.
9. The quantitation EMRTs are matched using the optimised parameters.

If a master identification file is used (`master` is set to `TRUE`, default is `FALSE`), the relevant actions that have already been executed when the file was created with [makeMaster](#) are not repeated here.

### Value

Invisibly returns an object of class `Synapter`. Used for its side effect of creating an html report of the run in `outputdir`.

### Author(s)

Laurent Gatto, Sebastian Gibb

### References

Bond N. J., Shliaha P.V., Lilley K.S. and Gatto L. (2013) J. Prot. Research.

### Examples

```
## Not run:
library(synapterdata)
data(synobj2)
output <- tempdir() ## a temporary directory
synergise2(object = synobj2, outputdir = output, outputfile = "synapter.html")
htmlReport <- paste0("file:/// ", file.path(output, "synapter.html")) ## the result report
browseURL(htmlReport) ## open the report with default browser

## End(Not run)
```

# Index

- \* **classes**
  - MasterFdrResults-class, [8](#)
  - MasterPeptides-class, [9](#)
  - Synapter, [13](#)
- \* **datasets**
  - synapterTinyData, [28](#)
- \* **internal**
  - Deprecated, [4](#)
- \* **methods**
  - plotFragmentMatching, [10](#)
- \* **package**
  - synapter-package, [2](#)
- allComb (MasterFdrResults-class), [8](#)
- allComb, MasterFdrResults-method (MasterFdrResults-class), [8](#)
- as.MSnSet.Synapter (Synapter), [13](#)
- bestComb (MasterFdrResults-class), [8](#)
- bestComb, MasterFdrResults-method (MasterFdrResults-class), [8](#)
- class:MasterFdrResults (MasterFdrResults-class), [8](#)
- class:MasterPeptides (MasterPeptides-class), [9](#)
- class:Synapter (Synapter), [13](#)
- createUniquePeptideDbRds, [3](#)
- Deprecated, [4](#)
- dim, Synapter-method (Synapter), [13](#)
- estimateMasterFdr, [4](#), [8](#)
- fileNames, MasterFdrResults-method (MasterFdrResults-class), [8](#)
- filterFragments (Synapter), [13](#)
- filterFragments, Synapter-method (Synapter), [13](#)
- filterIdentPepScore (Synapter), [13](#)
- filterIdentPepScore, Synapter-method (Synapter), [13](#)
- filterIdentPpmError (Synapter), [13](#)
- filterIdentPpmError, Synapter-method (Synapter), [13](#)
- filterIdentProtFpr (Synapter), [13](#)
- filterIdentProtFpr, Synapter-method (Synapter), [13](#)
- filterNonUniqueIdentMatches (Synapter), [13](#)
- filterNonUniqueIdentMatches, Synapter-method (Synapter), [13](#)
- filterNonUniqueMatches (Synapter), [13](#)
- filterNonUniqueMatches, Synapter-method (Synapter), [13](#)
- filterPeptideLength (Synapter), [13](#)
- filterPeptideLength, Synapter-method (Synapter), [13](#)
- filterQuantPepScore (Synapter), [13](#)
- filterQuantPepScore, Synapter-method (Synapter), [13](#)
- filterQuantPpmError (Synapter), [13](#)
- filterQuantPpmError, Synapter-method (Synapter), [13](#)
- filterQuantProtFpr (Synapter), [13](#)
- filterQuantProtFpr, Synapter-method (Synapter), [13](#)
- filterUniqueDbPeptides (Synapter), [13](#)
- filterUniqueDbPeptides, Synapter-method (Synapter), [13](#)
- filterUniqueIdentDbPeptides (Synapter), [13](#)
- filterUniqueIdentDbPeptides, Synapter-method (Synapter), [13](#)
- filterUniqueMatches (Synapter), [13](#)
- filterUniqueMatches, Synapter-method (Synapter), [13](#)
- filterUniqueQuantDbPeptides (Synapter), [13](#)

- filterUniqueQuantDbPeptides, Synapter-method (Synapter), [13](#)
- findEMRTs (Synapter), [13](#)
- findEMRTs, Synapter-method (Synapter), [13](#)
- fragmentMatching, [10](#)
- fragmentMatching (Synapter), [13](#)
- fragmentMatching, Synapter-method (Synapter), [13](#)
- fragmentMatchingPerformance (Synapter), [13](#)
- fragmentMatchingPerformance, Synapter-method (Synapter), [13](#)
- getBestGridParams (Synapter), [13](#)
- getBestGridParams, Synapter-method (Synapter), [13](#)
- getBestGridValue (Synapter), [13](#)
- getBestGridValue, Synapter-method (Synapter), [13](#)
- getEMRTtable (Synapter), [13](#)
- getEMRTtable, Synapter-method (Synapter), [13](#)
- getFragmentMatchingPpmTolerance (Synapter), [13](#)
- getFragmentMatchingPpmTolerance, Synapter-method (Synapter), [13](#)
- getGrid (Synapter), [13](#)
- getGrid, Synapter-method (Synapter), [13](#)
- getGridDetails (Synapter), [13](#)
- getGridDetails, Synapter-method (Synapter), [13](#)
- getIdentificationFragments (Synapter), [13](#)
- getIdentificationFragments, Synapter-method (Synapter), [13](#)
- getIdentPpmError (Synapter), [13](#)
- getIdentPpmError, Synapter-method (Synapter), [13](#)
- getImDiff (Synapter), [13](#)
- getImDiff, Synapter-method (Synapter), [13](#)
- getLog (Synapter), [13](#)
- getLog, Synapter-method (Synapter), [13](#)
- getLowessSpan (Synapter), [13](#)
- getLowessSpan, Synapter-method (Synapter), [13](#)
- getPepNumbers, [5](#)
- getPepNumbers (Synapter), [13](#)
- getPepNumbers, Synapter-method (Synapter), [13](#)
- getPepScoreFdr (Synapter), [13](#)
- getPepScoreFdr, Synapter-method (Synapter), [13](#)
- getPpmErrorQs (Synapter), [13](#)
- getPpmErrorQs, Synapter-method (Synapter), [13](#)
- getProtFpr (Synapter), [13](#)
- getProtFpr, Synapter-method (Synapter), [13](#)
- getQuantitationSpectra (Synapter), [13](#)
- getQuantitationSpectra, Synapter-method (Synapter), [13](#)
- getQuantPpmError (Synapter), [13](#)
- getQuantPpmError, Synapter-method (Synapter), [13](#)
- getRtNsd (Synapter), [13](#)
- getRtNsd, Synapter-method (Synapter), [13](#)
- getRtQs (Synapter), [13](#)
- getRtQs, Synapter-method (Synapter), [13](#)
- inputFiles (Synapter), [13](#)
- inputFiles, Synapter-method (Synapter), [13](#)
- inspectPeptideScores, [5](#)
- isCurrent, Synapter-method (Synapter), [13](#)
- loess, [15–17](#), [31](#), [34](#), [35](#)
- makeMaster, [5](#), [6](#), [9](#), [14](#), [29](#), [31](#), [32](#), [35](#)
- masterFdr (MasterFdrResults-class), [8](#)
- masterFdr, MasterFdrResults-method (MasterFdrResults-class), [8](#)
- MasterFdrResults, [5](#)
- MasterFdrResults (MasterFdrResults-class), [8](#)
- MasterFdrResults-class, [8](#)
- MasterPeptides, [7](#), [14](#), [29](#), [32](#)
- MasterPeptides-class, [9](#)
- mergePeptides (Synapter), [13](#)
- mergePeptides, Synapter-method (Synapter), [13](#)
- modelIntensity (Synapter), [13](#)
- modelIntensity, Synapter-method (Synapter), [13](#)
- modelRt (Synapter), [13](#)
- modelRt, Synapter-method (Synapter), [13](#)
- MSnExp, [7](#), [23](#)
- MSnSet, [11–13](#), [27](#), [28](#)
- par, [10](#)

- performance (Synapter), [13](#)
- performance, Synapter-method (Synapter), [13](#)
- performance2 (Synapter), [13](#)
- performance2, Synapter-method (Synapter), [13](#)
- plot, MasterFdrResults, missing-method (MasterFdrResults-class), [8](#)
- plot.default, [10](#)
- plot.qvalue, [22](#)
- plotCumulativeNumberOfFragments (Synapter), [13](#)
- plotCumulativeNumberOfFragments, Synapter-method (Synapter), [13](#)
- plotEMRTtable (Synapter), [13](#)
- plotEMRTtable, Synapter-method (Synapter), [13](#)
- plotFdr (Synapter), [13](#)
- plotFdr, Synapter-method (Synapter), [13](#)
- plotFeatures (Synapter), [13](#)
- plotFeatures, Synapter-method (Synapter), [13](#)
- plotFragmentMatching, [10](#), [22](#)
- plotFragmentMatching, Synapter-method (plotFragmentMatching), [10](#)
- plotFragmentMatchingPerformance (Synapter), [13](#)
- plotFragmentMatchingPerformance, Synapter-method (Synapter), [13](#)
- plotGrid (Synapter), [13](#)
- plotGrid, Synapter-method (Synapter), [13](#)
- plotIntensity (Synapter), [13](#)
- plotIntensity, Synapter-method (Synapter), [13](#)
- plotPepScores, [5](#)
- plotPepScores (Synapter), [13](#)
- plotPepScores, Synapter-method (Synapter), [13](#)
- plotPpmError (Synapter), [13](#)
- plotPpmError, Synapter-method (Synapter), [13](#)
- plotRt (Synapter), [13](#)
- plotRt, Synapter-method (Synapter), [13](#)
- plotRtDiffs (Synapter), [13](#)
- plotRtDiffs, Synapter-method (Synapter), [13](#)
- requantify (requantify, MSnSet-method), [11](#)
- requantify, MSnSet-method, [11](#)
- rescaleForTop3 (rescaleForTop3, MSnSet, MSnSet-method), [12](#)
- rescaleForTop3, MSnSet, MSnSet-method, [12](#)
- rescueEMRTs (Synapter), [13](#)
- rescueEMRTs, Synapter-method (Synapter), [13](#)
- searchGrid (Synapter), [13](#)
- searchGrid, Synapter-method (Synapter), [13](#)
- setBestGridParams (Synapter), [13](#)
- setBestGridParams, Synapter-method (Synapter), [13](#)
- setFragmentMatchingPpmTolerance (Synapter), [13](#)
- setFragmentMatchingPpmTolerance, Synapter-method (Synapter), [13](#)
- setIdentPpmError (Synapter), [13](#)
- setIdentPpmError, Synapter-method (Synapter), [13](#)
- setImDiff (Synapter), [13](#)
- setImDiff, Synapter-method (Synapter), [13](#)
- setLowessSpan (Synapter), [13](#)
- setLowessSpan, Synapter-method (Synapter), [13](#)
- setPepScoreFdr (Synapter), [13](#)
- setPepScoreFdr, Synapter-method (Synapter), [13](#)
- setPpmError (Synapter), [13](#)
- setPpmError, Synapter-method (Synapter), [13](#)
- setProtFpr (Synapter), [13](#)
- setProtFpr, Synapter-method (Synapter), [13](#)
- setQuantPpmError (Synapter), [13](#)
- setQuantPpmError, Synapter-method (Synapter), [13](#)
- setRtNsd (Synapter), [13](#)
- setRtNsd, Synapter-method (Synapter), [13](#)
- show, MasterFdrResults-method (MasterFdrResults-class), [8](#)
- show, MasterPeptides-method (MasterPeptides-class), [9](#)
- show, Synapter-method (Synapter), [13](#)
- showFdrStats (Synapter), [13](#)

showFdrStats, Synapter-method  
    (Synapter), [13](#)  
Synapter, [3](#), [7](#), [8](#), [10](#), [13](#), [28–34](#)  
synapter (synapter-package), [2](#)  
synapter-defunct (Deprecated), [4](#)  
synapter-deprecated (Deprecated), [4](#)  
synapter-package, [2](#)  
synapterGUI, [7](#)  
synapterGUI (Deprecated), [4](#)  
synapterGuide, [26](#)  
synapterPlgsAgreement  
    (synapterPlgsAgreement, MSnSet-method),  
    [27](#)  
synapterPlgsAgreement, MSnSet-method,  
    [27](#)  
synapterTiny (synapterTinyData), [28](#)  
synapterTinyData, [28](#)  
synergise, [3](#), [4](#), [7](#), [29](#)  
synergise1, [34](#)  
synergise1 (synergise), [29](#)  
synergise2, [32](#)  
synergize (synergise), [29](#)  
synergize1 (synergise), [29](#)  
synergize2 (synergise2), [32](#)  
  
updateObject, Synapter-method  
    (Synapter), [13](#)  
  
validObject, Synapter-method (Synapter),  
    [13](#)  
  
writeFragmentLibrary (makeMaster), [6](#)  
writeFragmentLibrary, MasterPeptides, character-method  
    (makeMaster), [6](#)  
writeIdentPeptides (Synapter), [13](#)  
writeIdentPeptides, Synapter-method  
    (Synapter), [13](#)  
writeMasterPeptides (makeMaster), [6](#)  
writeMasterPeptides, MasterPeptides, character-method  
    (makeMaster), [6](#)  
writeMatchedEMRTs, [27](#)  
writeMatchedEMRTs (Synapter), [13](#)  
writeMatchedEMRTs, Synapter-method  
    (Synapter), [13](#)  
writeMergedPeptides (Synapter), [13](#)  
writeMergedPeptides, Synapter-method  
    (Synapter), [13](#)  
writeQuantPeptides (Synapter), [13](#)  
writeQuantPeptides, Synapter-method  
    (Synapter), [13](#)